

# FMH RapidScreen

Results you can trust



**Trust is a precious thing.** Immucor's dedication to innovation means we are committed to delivering enhanced products that improve your laboratory's operations.

**FMH RapidScreen** offers a more simplified feto-maternal hemorrhage screening test to give you peace of mind in generating an improved result for the protection of mothers and their newborns.



Transfuse | Transplant | Transform **a life**



# FMH RapidScreen

Immucor's faster, easier feto-maternal hemorrhage screening test.

**FMH RapidScreen offers key improvements over currently available fetal screen kits**

## Faster, easier test procedure

- 5 minute, room temperature incubation
- Fewer microscopic fields simplifies the reading and interpretation of test results
- Up to 35% decrease in total test time compared to currently marketed tests
- Minimizes the need to carry out a quantitative procedure in all cases and can be performed within the blood bank laboratory

## More consistent performance

- Monoclonal IgM Anti-D antibodies replace chemically modified polyclonal antibodies for improved performance
- Elimination of enzyme-treatment to the indicator cells aids in reduction of false positive test results

## Ordering Information

Product Code	Description	Size	Use
0007893	FMH RapidScreen	1 Kit (30-35 tests)	Detection of D+ red cells in D- mothers post-delivery

To learn more or to set up a consultation with one of our blood bank business managers, call **855.IMMUCOR** (855.466.8267) or visit **[www.immucor.com](http://www.immucor.com)**.

### In Canada:

t 1.800.565.0653

e [DBLCustomerService@immucor.com](mailto:DBLCustomerService@immucor.com)

